



MedGyn

Innovative Products for Ob/Gyn Professionals

Model AL-106

Digital Video Colposcope and Imaging Data System

Model AL-106



Model AL-106SA



- High-resolution CCD and fast auto-focusing technology
- Maintenance-free and cold LED light source
- Real-time display and magnification

Model AL-106

Specifications

- Lens Resolution** ● 1/4" Super HAD color CCD, 630,000 pixels
- Horizontal Resolution** ● 470 lines
- Output Interface** ● S-Video + Video
- Working Distance** ● Normal 200-300mm (Autozoom)
- Magnification** ● x1 ~28x (typical)
- Field of View** ● 12~120mm
- Depth of View** ● 5~180mm
- Light Source** ● Multi-point LED light source
- Colors** ● White and green light, with bright or dark selection
- Image Freeze** ● Button Selection
- Focusing Method** ● Fast automatic and manual focusing
- Image Display** ● Color LCD monitor or Laptop
- Computer** ● Video image compatible
- Printer** ● Optional
- Continuous Acquisition Capacity** ● 1-30 frames
- Power Requirement** ● AC100-240V
- Working Temperature** ● 5°C-45°C
- Working Humidity** ● 30-80%RH

MedGyn Products Inc/International Inc
328 N. Eisenhower Lane
Lombard, IL U.S.A. 60148

www.medgyn.com email medgyn@medgyn.com
630.627.4105 toll free 800.451.9667
fax 630.627.0127



MedGyn

Innovative Products for Ob/Gyn Professionals