



MedGyn

Innovative Products for Ob/Gyn Professionals

Binocular Colposcope

Model AL-401
with 3 - Step
magnifications



AL 401 SERIES FEATURES

- Precision computer designed 10x optics
- Interchangeable eyepieces
- Individually adjustable eyepieces
- Electronically controlled illumination
- Smooth action controls
- Green filter

phone: 630.627.4105
toll-free US Sales: 800.451.9667
toll-free: 888.563.3496
fax: 630.627.0127

328 N. Eisenhower Lane
Lombard, IL U.S.A. 60148
email: medgyn@medgyn.com
web: <http://www.medgyn.com>


AL-401 Binocular Colposcope

General Specifications and Features

EASIER. TO OPERATE, INDIVIDUALLY MOVABLE BINOCULAR EYEPIECES
standard magnification included 10x

MedGyn AL-401 series

Colposcopes are the most versatile and highly developed. The standard features of the AL-401 are parfocal optics in 5x, 10x and 20x magnification, 300mm focal length, adjustable halogen light source, green filter, and interchangeable eyepieces with rubber cups.



MICROSCOPE		
Eye pieces	10X	
Magnifying Power	10X	
Total Magnification	5X, 10x, 20x	
Working Distance	300mm	
Field of Vision	20mm ø	with 10X eyepieces

STAND		
AL-401		
Tilting Adjustment	From Horizon	14° upward 45° downward

ILLUMINATION		
Lamp	12V, 20W	
Power Supply Unit	12V, 2.5A	
Filter	Built-in Green Filter for Internal Method	
Main Power	120VAC	
Illumination Field	50mm	

Net Weight	8kg	18lbs.
------------	-----	--------

328 N. Eisenhower Lane
Lombard, IL U.S.A. 60148

phone: 630.627.4105

fax: 630.627.0127

email: medgyn@medgyn.com

http://www.medgyn.com

toll-free US Sales: **800.451.9667**

888.563.3496



MedGyn

Innovative Products for Ob/Gyn Professionals